

KERN PES-A03

Loop for underfloor weighing

KERN



Category

Brand	KERN
Product category	Balance
Product group	Underfloor weighing

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	185×130×15 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Product information

GTIN/EAN number	4045761078260
REACH classification	no prohibited substances

Pictograms

STANDARD

